

Analytik + HPLC GmbH











DEM

HPLC capillary connectors with

1/32" (0.8 mm) OD or 1/16" (1.6 mm) OD made from Glass-PEEK (PEEK lined fused-silica-capillaries)

Glass T EETC (TEER III) ca Tasca silica capillaries/

Please request Technical Note 002e!





Please request Technical Note **003**e!









HPLC SYSTEM ACCESSORIES

Advantages:

- ideal replacement for stainless steel and PEEK capillaries
- suitable for capillary LC and LC-MS
- absolutely smooth inner surface (glass) and thus excellent flow characteristics
- no contamination
- no band broadening
- outstanding reproducibility
- no clogging
- metal-free
- suitable for high pressure applications
- connections readily made with the straight yet highly flexible elements

Large selection of capillary HPLC columns in a wide variety of dimensions and forms

HPLC columns with IDs from 50 μm to 800 μm, filled with any stationary phase.
The configuration of individual columns can be tailored to the customer's requirements

tailored to the customer's requirements, guaranteeing a perfect fit in form and function with all HPLC systems.

PEEK M4 fingertight fitting for Agilent Cap-LC and similar systems

Fingertight fitting with M4 threads for capillaries with 1/32" (0.8 mm) OD

Order Nr. 6941 5 pieces 65.00 Euro Order Nr. 6942 10 pieces 99.00 Euro

Guard column holders ID 20 and 30 mm

Fingertight as a result of new O-ring construction.

Special price: 39.00 Euro (each)

Replacement rotors for Rheodyne injectors Type 7125 and 7725 Order Nr. I 7125-047 Type 7010 Order Nr. I 7010-039